

Suppressor of cytokine signaling 3

Catalog No: #AP88182



Package Size: #AP88182-1 50ug #AP88182-2 100ug #AP88182-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Suppressor of cytokine signaling 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SOCS-3,Cytokine-inducible SH2 protein 3,CIS-3,Protein EF-10,Cis3,Cish3
Accession No.	Uniprot ID:O35718
Uniprot	O35718
GeneID	12702;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Socs3

Note: This product is for in vitro research use only