

SAM pointed domain-containing Ets transcription factor

Catalog No: #AP88186

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP88186-1 50ug #AP88186-2 100ug #AP88186-3 1mg

Description

Product Name	SAM pointed domain-containing Ets transcription factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Prostate epithelium-specific Ets transcription factor,Prostate-specific Ets,Prostate-derived Ets factor,Pdef,Pse
Accession No.	Uniprot ID:Q9WTP3
Uniprot	Q9WTP3
GeneID	30051;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Spdef

Note: This product is for in vitro research use only