

Angiopoietin-1 receptor

Catalog No: #AP88187



Package Size: #AP88187-1 50ug #AP88187-2 100ug #AP88187-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiopoietin-1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Endothelial tyrosine kinase, Tyrosine kinase with Ig and EGF homology domains-2, Tyrosine-protein kinase receptor TIE-2, CD202b, TIE-2, TIE2
Accession No.	Uniprot ID:Q06807
Uniprot	Q06807
GeneID	280939;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TEK

Note: This product is for in vitro research use only