

Tenascin

Catalog No: #AP88191



Package Size: #AP88191-1 50ug #AP88191-2 100ug #AP88191-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tenascin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TN,Cytotactin,GMEM,GP 150-225,Glioma-associated-extracellular matrix antigen,Hexabrachion,JI,Myotendinous antigen,Neuronectin,Tenascin-C,TN-C,HXB
Accession No.	Uniprot ID:P24821
Uniprot	P24821
GeneID	3371;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNC

Note: This product is for in vitro research use only