

TSC22 domain family protein 1

Catalog No: #AP88193



Package Size: #AP88193-1 50ug #AP88193-2 100ug #AP88193-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TSC22 domain family protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cerebral protein 2,Regulatory protein TSC-22,TGFB-stimulated clone 22 homolog,Transforming growth factor beta-1-induced transcript 4 protein,KIAA1994,TGFB114,TSC22,hucep-2
Accession No.	Uniprot ID:Q15714
Uniprot	Q15714
GeneID	8848;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TSC22D1

Note: This product is for in vitro research use only