

E3 ubiquitin-protein ligase SIAH2

Catalog No: #AP88198



Package Size: #AP88198-1 50ug #AP88198-2 100ug #AP88198-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase SIAH2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	RING-type E3 ubiquitin transferase SIAH2,Seven in absentia homolog 2,Siah-2
Accession No.	Uniprot ID:Q8R4T2
Uniprot	Q8R4T2
GeneID	140593;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Siah2

Note: This product is for in vitro research use only