

Solute carrier family 22 member 1

Catalog No: #AP88202



Package Size: #AP88202-1 50ug #AP88202-2 100ug #AP88202-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Solute carrier family 22 member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Organic cation transporter 1,hOCT1,OCT1
Accession No.	Uniprot ID:O15245
Uniprot	O15245
GeneID	6580;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC22A1

Note: This product is for in vitro research use only