

# Succinate-semialdehyde dehydrogenase, mitochondrial

Catalog No: #AP88208

Package Size: #AP88208-1 50ug #AP88208-2 100ug #AP88208-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Succinate-semialdehyde dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Aldehyde dehydrogenase family 5 member A1,NAD(+)-dependent succinic semialdehyde dehydrogenase, Ssadh
Accession No.	Uniprot ID:P51650
Uniprot	P51650
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Aldh5a1

Note: This product is for in vitro research use only