

Somatotropin

Catalog No: #AP88210



Package Size: #AP88210-1 50ug #AP88210-2 100ug #AP88210-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Somatotropin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Growth hormone,GH,Growth hormone 1,Pituitary growth hormone
Accession No.	Uniprot ID:P01241
Uniprot	P01241
GeneID	2688;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GH1

Note: This product is for in vitro research use only