

VPS10 domain-containing receptor SorCS1

Catalog No: #AP88211



Package Size: #AP88211-1 50ug #AP88211-2 100ug #AP88211-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VPS10 domain-containing receptor SorCS1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hSorCS,SORCS
Accession No.	Uniprot ID:Q8WY21
Uniprot	Q8WY21
GeneID	114815;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SORCS1

Note: This product is for in vitro research use only