

Stanniocalcin-2

Catalog No: #AP88212



Package Size: #AP88212-1 50ug #AP88212-2 100ug #AP88212-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Stanniocalcin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	STC-2
Accession No.	Uniprot ID:O88452
Uniprot	O88452
GeneID	20856;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Stc2

Note: This product is for in vitro research use only