

Pulmonary surfactant-associated protein D

Catalog No: #AP88214



Package Size: #AP88214-1 50ug #AP88214-2 100ug #AP88214-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pulmonary surfactant-associated protein D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PSP-D,CP4,Lung surfactant protein D,Sftp4
Accession No.	Uniprot ID:P35248
Uniprot	P35248
GeneID	25350;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sftpd

Note: This product is for in vitro research use only