

Signal recognition particle 54 kDa protein

Catalog No: #AP88217



Package Size: #AP88217-1 50ug #AP88217-2 100ug #AP88217-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Signal recognition particle 54 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SRP54
Accession No.	Uniprot ID:Q2T9U1
Uniprot	Q2T9U1
GeneID	614017;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SRP54

Note: This product is for in vitro research use only