

Tumor necrosis factor receptor superfamily member 11A

Catalog No: #AP88218

Package Size: #AP88218-1 50ug #AP88218-2 100ug #AP88218-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 11A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Osteoclast differentiation factor receptor,ODFR,Receptor activator of NF-KB,CD265,RANK
Accession No.	Uniprot ID:Q9Y6Q6
Uniprot	Q9Y6Q6
GeneID	8792;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFRSF11A

Note: This product is for in vitro research use only