

Tumor necrosis factor receptor superfamily member 16

Catalog No: #AP88219

Package Size: #AP88219-1 50ug #AP88219-2 100ug #AP88219-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Gp80-LNGFR,Low affinity neurotrophin receptor p75NTR,Low-affinity nerve growth factor receptor,NGF receptor,p75 ICD,TNFRSF16
Accession No.	Uniprot ID:P18519
Uniprot	P18519
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NGFR

Note: This product is for in vitro research use only