

Tumor necrosis factor receptor superfamily member 10D

Catalog No: #AP88221

Package Size: #AP88221-1 50ug #AP88221-2 100ug #AP88221-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 10D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Decoy receptor 2,DcR2,TNF-related apoptosis-inducing ligand receptor 4,TRAIL receptor 4,TRAIL receptor with a truncated death domain,CD264,DCR2,TRAILR4,TRUNDD,UNQ251/PRO288
Accession No.	Uniprot ID:Q9UBN6
Uniprot	Q9UBN6
GeneID	8793;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFRSF10D

Note: This product is for in vitro research use only