

Tripartite motif-containing protein 72

Catalog No: #AP88224



Package Size: #AP88224-1 50ug #AP88224-2 100ug #AP88224-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tripartite motif-containing protein 72
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Mitsugumin-53,Mg53,Mg53
Accession No.	Uniprot ID:A0JPQ4
Uniprot	A0JPQ4
GeneID	365377;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Trim72

Note: This product is for in vitro research use only