

Tribbles homolog 3

Catalog No: #AP88225



Package Size: #AP88225-1 50ug #AP88225-2 100ug #AP88225-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tribbles homolog 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TRB-3, Neuronal cell death-inducible putative kinase, SINK, p65-interacting inhibitor of NF-kappa-B, C20orf97, NIPK, SKIP3, TRB3
Accession No.	Uniprot ID: Q96RU7
Uniprot	Q96RU7
GeneID	57761;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TRIB3

Note: This product is for in vitro research use only