

Starch-binding domain-containing protein 1

Catalog No: #AP88228



Package Size: #AP88228-1 50ug #AP88228-2 100ug #AP88228-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Starch-binding domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Genethonin-1,Glycophagy cargo receptor stbd1,D5Erttd593e
Accession No.	Uniprot ID:Q8C7E7
Uniprot	Q8C7E7
GeneID	52331;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Stbd1

Note: This product is for in vitro research use only