

Semaphorin-3E

Catalog No: #AP88230



Package Size: #AP88230-1 50ug #AP88230-2 100ug #AP88230-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Semaphorin-3E
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	KIAA0331,SEMAH
Accession No.	Uniprot ID:O15041
Uniprot	O15041
GeneID	9723;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SEMA3E

Note: This product is for in vitro research use only