

Selenoprotein P

Catalog No: #AP88231



Package Size: #AP88231-1 50ug #AP88231-2 100ug #AP88231-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Selenoprotein P
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SeP, Plasma selenoprotein P, Selp, Sepp1
Accession No.	Uniprot ID:P70274
Uniprot	P70274
GeneID	20363;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Selenop

Note: This product is for in vitro research use only