

Pulmonary surfactant-associated protein A

Catalog No: #AP88233



Package Size: #AP88233-1 50ug #AP88233-2 100ug #AP88233-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pulmonary surfactant-associated protein A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PSAP,Sftp-1,Sftp1,Sftpa
Accession No.	Uniprot ID:P08427
Uniprot	P08427
GeneID	24773;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sftpa1

Note: This product is for in vitro research use only