

Serotonin N-acetyltransferase

Catalog No: #AP88234



Package Size: #AP88234-1 50ug #AP88234-2 100ug #AP88234-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serotonin N-acetyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Serotonin acetylase,Aralkylamine N-acetyltransferase,AA-NAT,SNAT
Accession No.	Uniprot ID:P79774
Uniprot	P79774
GeneID	396066;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AANAT

Note: This product is for in vitro research use only