

Selenium-binding protein 1

Catalog No: #AP88237



Package Size: #AP88237-1 50ug #AP88237-2 100ug #AP88237-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Selenium-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	56 kDa selenium-binding protein, SBP56, SBP
Accession No.	Uniprot ID:Q2KJ32
Uniprot	Q2KJ32
GeneID	510154;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SELENBP1

Note: This product is for in vitro research use only