

Acetyl-CoA acetyltransferase, mitochondrial

Catalog No: #AP88241



Package Size: #AP88241-1 50ug #AP88241-2 100ug #AP88241-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Acetyl-CoA acetyltransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acetoacetyl-CoA thiolase,T2,ACAT,MAT
Accession No.	Uniprot ID:P24752
Uniprot	P24752
GeneID	38;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ACAT1

Note: This product is for in vitro research use only