

Neurokinin-B

Catalog No: #AP88242



Package Size: #AP88242-1 50ug #AP88242-2 100ug #AP88242-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurokinin-B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	NKB, Neuromedin-K, NKNB
Accession No.	Uniprot ID:P67934
Uniprot	P67934
GeneID	492314;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TAC3

Note: This product is for in vitro research use only