

## Ras-related protein M-Ras

Catalog No: #AP88243



Package Size: #AP88243-1 50ug #AP88243-2 100ug #AP88243-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ras-related protein M-Ras
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ras-related protein R-Ras3,RRAS3
Accession No.	Uniprot ID:O14807
Uniprot	O14807
GeneID	22808;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MRAS

Note: This product is for in vitro research use only