

## Thyroid hormone receptor alpha

Catalog No: #AP88244



Package Size: #AP88244-1 50ug #AP88244-2 100ug #AP88244-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Thyroid hormone receptor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Nuclear receptor subfamily 1 group A member 1,c-erbA-1,c-erbA-alpha,C-erba-alpha,Nr1a1
Accession No.	Uniprot ID:P63058
Uniprot	P63058
GeneID	21833;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Thra

Note: This product is for in vitro research use only