

Titin

Catalog No: #AP88245



Package Size: #AP88245-1 50ug #AP88245-2 100ug #AP88245-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Titin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Connectin,Rhabdomyosarcoma antigen MU-RMS-40.14
Accession No.	Uniprot ID:Q8WZ42
Uniprot	Q8WZ42
GeneID	7273;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TTN

Note: This product is for in vitro research use only