

Serine/threonine-protein kinase SIK1

Catalog No: #AP88248

Package Size: #AP88248-1 50ug #AP88248-2 100ug #AP88248-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase SIK1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HRT-20,Myocardial SNF1-like kinase,Salt-inducible kinase 1,SIK-1,Serine/threonine-protein kinase SNF1-like kinase 1,Serine/threonine-protein kinase SNF1LK,Msk,Sik,Snf1lk
Accession No.	Uniprot ID:Q60670
Uniprot	Q60670
GeneID	17691;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sik1

Note: This product is for in vitro research use only