

Translationally-controlled tumor protein

Catalog No: #AP88249



Package Size: #AP88249-1 50ug #AP88249-2 100ug #AP88249-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Translationally-controlled tumor protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TCTP, Fortilin, Histamine-releasing factor, HRF, p23
Accession No.	Uniprot ID:P13693
Uniprot	P13693
GeneID	7178;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TPT1

Note: This product is for in vitro research use only