

Signal recognition particle 19 kDa protein

Catalog No: #AP88250



Package Size: #AP88250-1 50ug #AP88250-2 100ug #AP88250-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Signal recognition particle 19 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SRP19
Accession No.	Uniprot ID:Q3ZBG7
Uniprot	Q3ZBG7
GeneID	514960;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SRP19

Note: This product is for in vitro research use only