

Estrogen sulfotransferase, isoform 2

Catalog No: #AP88251

Package Size: #AP88251-1 50ug #AP88251-2 100ug #AP88251-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Estrogen sulfotransferase, isoform 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	EST-2,Estrone sulfotransferase,Sulfotransferase, estrogen-preferring,Ste
Accession No.	Uniprot ID:P52845
Uniprot	P52845
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ste2

Note: This product is for in vitro research use only