

Transcription factor 7-like 2

Catalog No: #AP88254



Package Size: #AP88254-1 50ug #AP88254-2 100ug #AP88254-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor 7-like 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HMG box transcription factor 4, T-cell-specific transcription factor 4, T-cell factor 4, Tcf4
Accession No.	Uniprot ID:Q924A0
Uniprot	Q924A0
GeneID	21416;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tcf712

Note: This product is for in vitro research use only