

NAD-dependent deacetylase sirtuin-3

Catalog No: #AP88256

Package Size: #AP88256-1 50ug #AP88256-2 100ug #AP88256-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD-dependent deacetylase sirtuin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SIR2-like protein 3,mSIR2L3,Sir2I3
Accession No.	Uniprot ID:Q8R104
Uniprot	Q8R104
GeneID	64384;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sirt3

Note: This product is for in vitro research use only