

Mothers against decapentaplegic homolog 3

Catalog No: #AP88258



Package Size: #AP88258-1 50ug #AP88258-2 100ug #AP88258-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mothers against decapentaplegic homolog 3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | MAD homolog 3, SMAD family member 3, SMAD 3, Madh3 |
| Accession No. | Uniprot ID:P84025 |
| Uniprot | P84025 |
| GeneID | 25631; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Smad3

Note: This product is for in vitro research use only