

Mothers against decapentaplegic homolog 5

Catalog No: #AP88259



Package Size: #AP88259-1 50ug #AP88259-2 100ug #AP88259-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAD homolog 5,Dwarfin-C,Dwf-C,SMAD family member 5,SMAD 5,Madh5,Msmad5
Accession No.	Uniprot ID:P97454
Uniprot	P97454
GeneID	17129;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Smad5

Note: This product is for in vitro research use only