

Thy-1 membrane glycoprotein

Catalog No: #AP88265



Package Size: #AP88265-1 50ug #AP88265-2 100ug #AP88265-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thy-1 membrane glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Thy-1 antigen,CD90,Thy-1
Accession No.	Uniprot ID:P01831
Uniprot	P01831
GeneID	21838;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Thy1

Note: This product is for in vitro research use only