

Signal transducer and activator of transcription 1-alpha/beta

Catalog No: #AP88266

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88266-1 50ug #AP88266-2 100ug #AP88266-3 1mg

Description

Product Name	Signal transducer and activator of transcription 1-alpha/beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Transcription factor ISGF-3 components p91/p84
Accession No.	Uniprot ID:P42224
Uniprot	P42224
GeneID	6772;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:STAT1

Note: This product is for in vitro research use only