

Trimethyllysine dioxygenase, mitochondrial

Catalog No: #AP88274



Package Size: #AP88274-1 50ug #AP88274-2 100ug #AP88274-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Trimethyllysine dioxygenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Epsilon-trimethyllysine 2-oxoglutarate dioxygenase, TML hydroxylase, TML-alpha-ketoglutarate dioxygenase, TML dioxygenase, Tmlh
Accession No.	Uniprot ID:Q91ZW6
Uniprot	Q91ZW6
GeneID	170898;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tmlhe

Note: This product is for in vitro research use only