

Sentan

Catalog No: #AP88275



Package Size: #AP88275-1 50ug #AP88275-2 100ug #AP88275-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sentan
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein S100-A1-like,S100 calcium-binding protein A1-like,S100a1l
Accession No.	Uniprot ID:Q8C9X1
Uniprot	Q8C9X1
GeneID	218739;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sntn

Note: This product is for in vitro research use only