

Sodium-dependent dopamine transporter

Catalog No: #AP88277



Package Size: #AP88277-1 50ug #AP88277-2 100ug #AP88277-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium-dependent dopamine transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	DA transporter,Solute carrier family 6 member 3
Accession No.	Uniprot ID:P27922
Uniprot	P27922
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC6A3

Note: This product is for in vitro research use only