

Tumor necrosis factor ligand superfamily member 4

Catalog No: #AP88283



Package Size: #AP88283-1 50ug #AP88283-2 100ug #AP88283-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	OX40 ligand, OX40L, CD252, Ox40I, Txgp1I
Accession No.	Uniprot ID:P43488
Uniprot	P43488
GeneID	22164;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tnfsf4

Note: This product is for in vitro research use only