

Tenascin-R

Catalog No: #AP88285



Package Size: #AP88285-1 50ug #AP88285-2 100ug #AP88285-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Tenascin-R |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | TN-R,Janusin,Neural recognition molecule J1-160/180,Restrictin |
| Accession No. | Uniprot ID:Q05546 |
| Uniprot | Q05546 |
| GeneID | 25567; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnr

Note: This product is for in vitro research use only