

TLR4 interactor with leucine rich repeats

Catalog No: #AP88288



Package Size: #AP88288-1 50ug #AP88288-2 100ug #AP88288-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TLR4 interactor with leucine rich repeats
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Leucine-rich repeat-containing protein KIAA0644,Kiaa0644
Accession No.	Uniprot ID:Q9DBY4
Uniprot	Q9DBY4
GeneID	66873;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tril

Note: This product is for in vitro research use only