

## Thioredoxin reductase 1, cytoplasmic

Catalog No: #AP88289



Package Size: #AP88289-1 50ug #AP88289-2 100ug #AP88289-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Thioredoxin reductase 1, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TR,Thioredoxin reductase TR1
Accession No.	Uniprot ID:O62768
Uniprot	O62768
GeneID	282388;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TXNRD1

Note: This product is for in vitro research use only