

## Mothers against decapentaplegic homolog 1

Catalog No: #AP88293



Package Size: #AP88293-1 50ug #AP88293-2 100ug #AP88293-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mothers against decapentaplegic homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MAD homolog 1, SMAD family member 1, SMAD 1
Accession No.	Uniprot ID:Q1JQA2
Uniprot	Q1JQA2
GeneID	540488;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: SMAD1

Note: This product is for in vitro research use only