

Tryptophan 2,3-dioxygenase

Catalog No: #AP88294



Package Size: #AP88294-1 50ug #AP88294-2 100ug #AP88294-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tryptophan 2,3-dioxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TDO, Tryptamin 2,3-dioxygenase, Tryptophan oxygenase, TO, Tryptophan pyrrolase, Tryptophanase
Accession No.	Uniprot ID:Q2KIQ5
Uniprot	Q2KIQ5
GeneID	530397;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TDO2

Note: This product is for in vitro research use only