

Suppressor of cytokine signaling 3

Catalog No: #AP88297



Package Size: #AP88297-1 50ug #AP88297-2 100ug #AP88297-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Suppressor of cytokine signaling 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SOCS-3
Accession No.	Uniprot ID:Q9BEG9
Uniprot	Q9BEG9
GeneID	282081;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SOCS3

Note: This product is for in vitro research use only