

Sushi-repeat-containing protein SRPX

Catalog No: #AP88303



Package Size: #AP88303-1 50ug #AP88303-2 100ug #AP88303-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sushi-repeat-containing protein SRPX
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DRS protein
Accession No.	Uniprot ID:Q9R0M3
Uniprot	Q9R0M3
GeneID	51795;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SrpX

Note: This product is for in vitro research use only